



# The PhysioLITE III



The **PhysioLITE III** is a compact, clinical grade low intensity shockwave therapy system that induces an intracellular and extracellular reactions to assist in the healing and recovery of soft tissue.

Shockwave technology theory works by a deep penetrating focused wave that activate cells, driving the formation of new cells that replace damaged tissues. Shockwave therapy accelerates blood flow to the treated area and may induce vasodilation and angiogenesis restoring and improving musculoskeletal tissue.

- Clinical Strength
- Deep Penetrating Focused Shockwave
- Adjustable Frequency
- Pain Free 5 Minute Treatment
- Rigid High-Quality Construction
- Easy to Administer Therapy

The **PhysioLITE III** is a therapeutic shockwave device that generates deep penetrating impulses into soft tissue. The pulse waves travel deep into the tissue and induce vasodilation, which helps increase blood flow to the treated area and is found to help relieve pain and accelerate healing. The **PhysioLITE III** is ergonomic, light weight, and equipped with a simplified user interface for easily administered treatments. There are 3 frequency settings (10Hz,15Hz and 20Hz) and adjustable time on the hand control. The **PhysioLITE III** offers multiple output intensities to accommodate multitude of therapy protocols. This device can provide a wide range of pain management treatments. The **PhysioLITE III** is supplied with, ultrasound gel, AC/DC adapter, instruction manual, quick start, guide, and a book of recommended therapy protocols.



### The PhysioLITE II FEATURES:

- Multiple frequency settings for mild or acute pain management
- Focused deep penetrating shockwaves
- Therapy treatment guide included with over 25 therapy protocols
- Built-in timer
- Bright intuitive control and display
- Rigid high-quality construction
- Clinical grade power output
- Selectable Hz
- Small, ergonomic design
- Easy to administer therapy
- Comes complete with device, ultrasound gel, AC/DC adapter, quick start guide and instruction manual and therapy protocol.
- One year warranty

### The PhysioLITE III TECHNICAL SPECIFICATIONS:

<b>Acoustic Frequency</b>	10 Hz, 15 Hz, 20 Hz
<b>Shockwave</b>	Continuous Focused Acoustic Wave
<b>Waveform</b>	Continuous / Pulsed
<b>Power Output</b>	35 mJ / 1 bar
<b>Treatment Time</b>	1-5 minutes